



| | |
|---------------------------------------|---|
| fitted with | LED - Energy efficiency category A+ |
| Illuminancefield @ 100 cm | Ø 18 cm |
| illuminance @ 100 cm | 100 000 lx |
| color rendering index (CRI) | > 95 |
| colour temperatur | 3700/4300/4700 K |
| work equipment | electronic converter |
| mains lead | direct power supply |
| luminaire body | polyamide (PA) |
| weight (net) | 15.9 kg |
| Power consumption | 40 W |
| usage | keypad at luminaire head |
| balance of articulation joints | spring |
| class of protection | I |
| Accessory | sterilizable handle, adapter for suspended ceiling |
| norms | EN 60601-1, EN 60601-2-41 |
| special features | 3 color shades, 4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free |
| Fastening | ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm, intermediate ceiling adapter (on request) |
| System of protection | IP 20 (43) |
| order number | D15885000 - pure white |



light distribution

