

VISIANO 20-2



reddot award 2014
winner

Derungs

Tested safety

medical device IEC 60601,
unique in this class



Quality arm system

durable, robust and secure application



Unique diamond optics

top colour homogeneity of light



Colour change 3 500 K / 4 500 K

optimal contrast sensitivity



high-quality components

more durable, cost-effective application



Spring-balanced friction joint

easy positioning and locking



Optimal light

Ø 21 cm, 60 000 lx / 0.5 m, Ra > 95, R13 > 93
4 dimming levels



360° turnover-axis

maximum availability of light

universal fastening
on ceiling, wall or rollerstand



- examination room
- maternity ward
- emergency room
- intensive care unit
- neonatology
- gynaecology
- general practitioner
- dermatologist
- otorhinolaryngology