

Spring Arm System

TENO 30 C D15983000

Waldmann **W**

ENGINEER OF LIGHT.



Derungs **D**
MEDICAL LIGHTING

Payload	Max. 30Nm
Cable diameter	Max. 8mm
Vertical movement angle	+45° up / -65° down
Length extension arm	800mm*
Length spring arm	720mm*
Articulated joint	Spring balanced
Power supply	Integrated in ceiling plate 100-240 V; 50/60 Hz, 40 Watt
Mains lead	Direct power supply
Weight (netto)	ca. 13.5 kg
Accessory	Intermediate ceiling adapter
Normen	EN 60601-1
Special features	Closed construction, very good hygienic properties
Fastening	Ceiling fixture, ceiling tube can be cut to length (125-120 cm)
Color	Pure white, other colors on request

*other lengths on request

