

# Spring Arm System

TENO 30 W D15984000

# Waldmann **W**

ENGINEER OF LIGHT.



<b>Payload</b>	Max. 30Nm
<b>Cable diameter</b>	Max. 8mm
<b>Vertical movement angle</b>	+45° up / -65° down
<b>Length extension arm</b>	800mm*
<b>Length spring arm</b>	720mm*
<b>Articulated joint</b>	Spring balanced
<b>Power supply</b>	Integrated in ceiling plate 100-240 V; 50/60 Hz, 40 Watt
<b>Mains lead</b>	Direct power supply
<b>Weight (netto)</b>	ca. 13 kg
<b>Normen</b>	EN 60601-1
<b>Special features</b>	Closed construction, very good hygienic properties
<b>Fastening</b>	Wall fixture
<b>Color</b>	Pure white, other colors on request

**Derungs** **D**  
MEDICAL LIGHTING

\*other lengths on request

