

Spring Arm System

TENO 30 F D15985000

Waldmann **W**

ENGINEER OF LIGHT.

Derungs **D**
MEDICAL LIGHTING



Payload	Max. 30Nm
Cable diameter	Max. 8mm
Vertical movement angle	+45° up / -65° down
Length extension arm	800mm*
Length spring arm	720mm*
Articulated joint	Spring balanced
Power supply	Integrated in ceiling plate 100-240 V; 50/60 Hz, 40 Watt
Mains lead	Direct power supply
Weight (netto)	ca. 17.5 kg
Normen	EN 60601-1
Special features	Closed construction, very good hygienic properties
Fastening	Wall fixture
Color	Pure white, other colors on request

*other lengths on request

