



Payload	Max. 10Nm
Cable diameter	Max. 8mm
Vertical movement angle	93°
Length lower arm	600mm*
Length upper arm	480mm*
Articulated joint	Spring and friction balanced
Power supply	Integrated in ceiling plate 100-240 V; 50/60 Hz, 20 Watt**
Mains lead	1.8m cable
Weight (netto)	ca. 1.4 kg
Accessory	Intermediate ceiling adapter
Norms	EN 60601-1
Special features	Closed construction, very good hygienic properties
Fastening	Ceiling mount, ceiling tube can be cut on site, (length: 40-116cm)
Color	Pure white, other colors on request

*other lengths on request

**alternative power supply on request

Derungs 
MEDICAL LIGHTING

