



**Derungs**   
MEDICAL LIGHTING

<b>Payload</b>	Max. 10Nm
<b>Cable diameter</b>	Max. 8mm
<b>Vertical movement angle</b>	93°
<b>Length lower arm</b>	600mm*
<b>Length upper arm</b>	480mm*
<b>Articulated joint</b>	Spring and friction balanced
<b>Power supply</b>	Integrated in lower arm 100-240 V; 50/60 Hz, 20 Watt**
<b>Mains lead</b>	1.8m cable
<b>Weight (netto)</b>	ca. 1.4 kg
<b>Accessory</b>	Roller stand, rail mount, table clamp, universal mount, Wall mount
<b>Norms</b>	EN 60601-1
<b>Special features</b>	Closed construction, very good hygienic properties
<b>Fastening</b>	Fastening Pin (Ø16mm)
<b>Color</b>	Pure white, other colors on request

\*other lengths on request

\*\*alternative power supply on request

